INDIRECT FLUORESCENCE ASSAY FOR IgG ANTINUCLEAR ANTIBODY HEp2

For In Vitro Diagnostic Use

I-ANA01G 120 Tests

Intended Use

The SCIMEDX Indirect Fluorescence Assay (IFA) for IgG Antinuclear Antibody (ANA) HEp2 is intended for the qualitative and semi-quantitative detection of IgG (Immunoglobulin G) antinuclear antibodies in human serum. Detection of ANA in humans can be used as an aid in the diagnosis of autoimmune connective tissue diseases.

Introduction and Summary of Test Procedures

Systemic autoimmune diseases are characterized by a variety of autoantibodies specific for nuclear and cytoplasmic molecules involved in DNA replication, DNA transcription and messenger RNA (mRNA) translation. The autoantibodies are produced as a result of an imbalance in the natural immunoregulatory network, most likely caused by infectious agents (1-3). The manifestations of systemic autoimmunity are either due to the direct effects of these autoantibodies or to antigen-antibody deposition (1).

Circulating serum autoantibodies to nuclear antigens (ANA) are present in a number of diseases such as systemic lupus erythematosus (SLE), rheumatoid arthritis, discoid lupus, scleroderma, Sjögren's syndrome, liver disease, ulcerative colitis and myasthenia gravis. Approximately 5% of people without disease symptoms have significant levels of ANA (4, 5). The indirect fluorescent antibody (IFA) technique of Coons (6) was adapted to ANA testing by Friou (7, 8). It has since become the method of choice for screening and quantitation of ANA (9). The use of a tissue culture substrate, HEp2, provides a clear spectrum of nuclear staining patterns, unavailable with other substrates that may be associated with specific disease states (10-15). In addition, the HEp2 cell line is the recommended substrate for detection of anti-centromere antibody which is highly selective for the CREST variant of progressive systemic sclerosis (10).

Mixed patterns of ANA may commonly be observed with rheumatic disease sera, reflecting the simultaneous occurrence of multiple ANAs in the same serum (9). The IFA technique is very effective in localizing antigens in the nucleus, nucleoplasm or cytoplasm (9).

Principle of the Test

SCIMEDX fluorescent antibody assays use the indirect method of antibody detection and titer determination. Diluted patient serum samples are applied to cultured cells containing inactivated antigens provided on paint delineated wells on glass microscope slides. During a 30 minute incubation, antibody specific for nuclear antigens forms an antigen/antibody complex with the nuclear antigens in the cells. In a brief washing step, nonspecific antibody and other unreacted serum proteins are eliminated. Fluorescein-conjugated goat anti-human IgG is then applied to the wells of the glass slide. The anti-IgG conjugate combines with human IgG, if present, during a 30 minute incubation. After a brief wash to remove unreacted conjugate, the slides are viewed by fluorescence microscopy. A positive antibody reaction is denoted by bright green fluorescence at the antigen sites.

Materials Furnished and Storage Conditions

ANA HEp2 Antigen Slides: Slides of HEp2 cells on each glass well. The slides are ready for use after removal from protective pouch. Store at 2-8°C. Slides are stable at this storage condition until the expiration date stated on the pouch label.

ANA IgG Positive Control: Each vial contains 0.5 ml ANA IgG positive human control. This component is a ready for use liquid. Store at 2-8°C. Liquid positive control is stable at this storage condition until the expiration date stated on the vial label.

ANA IgG Negative Control: Each vial contains 0.5 ml ANA IgG negative human control. This component is a ready for use liquid. Store at 2-8°C. Liquid negative control is stable at this storage condition until the expiration date stated on the vial label.

Fluorescein Conjugate: Each vial contains 1.5 ml fluorescein conjugated goat (inactivated) antihuman IgG (heavy and light chain) with Evans Blue and Rhodamine counterstains. The fluorescein conjugate is a conjugation of affinity chromatography purified antihuman IgG with fluorescein isothiocyanate (FITC). Adding Evans Blue and Rhodamine counterstains to the conjugate masks nonspecific fluorescence of the tissue culture cells. This component is a ready for use liquid. Store at 2-8°C. Liquid conjugate is stable at this storage condition until the expiration date stated on the vial label.

Coverslip Mounting Media: Each vial contains 2.0 ml phosphate buffered glycerol with fade retardant. This component is a ready for use liquid. Storage temperature may range from refrigerated to room temperatures (2-30°C). Mounting media is stable at either storage condition until the expiration date stated on the vial label.

Phosphate Buffered Saline (PBS): Each aluminized sealed packet of powdered buffer makes one liter of 1X PBS. Storage temperatures may range from refrigerated to room temperature (2-30°C). Add the entire contents of a PBS packet to one liter of freshly prepared distilled or deionized water. Note: Addition of the salts while rapidly stirring the water will facilitate solubilization. Store PBS as a solution at 2-8°C.

Special Blotters: Absorbent blotters have pre-cut holes for use in drying the slide mask. Storage temperatures may range from refrigerated to room temperature (2-30°C).

Additional Materials Required

- Test tubes, racks, pipettes, microtiter plates and safety pipetting devices for making sample dilutions
- 37°C incubator
- Moist chamber for incubating slides
- Slide holder rack and staining dish for washing slides
- Coverslips: 22 X 50 mm No. 1 thickness glass

Fluorescence microscope: A fluorescence microscope equipped with the following was used to calibrate the controls and conjugate:

- 10X eyepiece
- 16X or 40X objectives
- Epi-illuminator with 50W halogen lamp
- FITC-excitation filter KP490
- Yellow absorbing filter K530
- Red suppression filter BG38

The fluorescein label has an excitation peak of 490 nm and an emission peak of 520 nm. Differences in endpoint reactivities and fluorescence intensities may be due to the type and condition of the fluorescence equipment used in your laboratory.

General Precautions

No US Standard of Potency. For *in vitro* diagnostic use only.

- Store all kit components at their recommended or suggested temperatures. **Do not freeze.**
- Do not use the components beyond the stated expiration date of each component.
- SCIMEDX optimizes all of the active components in each lot of its IFA kits as a unit. Do not mix components from different lots or from different sources.
- The controls and conjugate contain 0.095% sodium azide that, if allowed to accumulate, can form explosive compounds in lead and/or copper plumbing. Flush drains thoroughly if used to dispose of these materials.
- Use separate pipette tips for each sample and reagent to avoid cross contamination.
- Reagents should be inspected for evidence of bacterial or fungal contamination.

Antigen slides: All IFA antigen slides have fixed cells that do not contain any viable infectious agents. However, good laboratory practices (GLP) require careful slide handling and disposal as with any other potentially biohazardous laboratory material.

- Do not remove the slides from their protective pouch until ready for use.
- The reported concentration of anti-ANA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.
- Do not reuse substrate slide.

Human Controls: The human controls in these kits have all been tested for hepatitis B surface antigen (HBsAg) and human immunodeficiency virus (HIV) antibody by FDA licensed methods and found to be non-reactive. However, no test system can ensure the absence of these agents. Handle all human serum components, including those received in your laboratory for testing, as potentially biohazardous. For guidelines refer to the Centers for Disease Control/National Institutes of Health Manual, "Biosafety in Microbiological and Biomedical Laboratories," 1993.

Xn-Harmful Substance US

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		<ul style="list-style-type: none"> • Wash thoroughly after handling. • Wear eye/face protection
Signal Word	WARNING	<ul style="list-style-type: none"> • IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. • If eye irritation persists: Get medical advice/attention.
Hazard Statement	Causes serious eye irritation.	<ul style="list-style-type: none"> • Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

EU

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.
Signal Word	WARNING	P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.
Hazard Statement	H319 Causes serious eye irritation.	P337+P313 If irritation persists, get medical advice/attention.

Control and Conjugate Safety Precaution:

Concentration of sodium azide in these components is classified as Harmful and subject to the following risk phrase: "Harmful if swallowed and contact with acids liberates very toxic gas."

Specimen Collection, Storage and Limitations of Test Samples

- Separate aseptically collected serum from the red blood cells and store frozen (-10°C or lower) until ready for testing. Avoid repeated freezing and thawing. For recommendations for collecting and storing blood specimens refer to the NCCLS publication, "Approved Standard – Procedures for the Handling and Processing of Blood Specimens, H18A.
- If desired, store fresh liquid serum samples at 2° to 8°C for up to one week without loss of antibody activity.
- Do not use hemolyzed, lipemic or icteric samples.
- Do not use visibly contaminated samples.
- Serum samples must not be heat inactivated prior to use.
- Do not add anticoagulants or preservatives to serum samples.

IFA Procedure

1. For qualitative IgG antibody determination, prepare a 1:40 screening dilution of each test sample in PBS. Prepare all dilutions in a minimum volume of 0.10 ml with PBS as the diluent.
2. For semi-quantitative titration of sera, prepare two-fold serial dilutions of the serum sample in PBS, starting with a 1:40 dilution, and adding equal volumes of diluted serum or plasma and PBS for each consecutive dilution.
3. Remove slides from protective pouch and apply 1 drop (approximately 20 µl) of the diluted test sample(s) to each well. Add sufficient volume to completely cover each well, but cross-mixing of contents between wells should not occur. Note: Each day's test run requires one well each for positive control, negative control, and PBS (conjugate control). Incubate the slides in a moist chamber for 30 minutes at 37°C.
4. Rinse the slides in a light stream of buffer. Avoid directing the stream at the wells.
5. Wash the slides for 10 minutes with a change of PBS solution after 5 minutes. Agitate the slides by moving the rack up and down in the buffer.
6. Blot the paint mask surrounding the test wells with the special blotters.
7. Apply one drop of the ready to use conjugate to each test well.
8. Repeat steps 4 (incubation), 5 (PBS rinse), 6 (10 minute PBS wash), and 7 (blot).
9. Apply the glycerol mounting media and 22 X 50 mm glass coverslip.
10. Observe the reactivity under fluorescence microscopy using 20-40X magnification. For best results, examine slides immediately after completion of the test. To obtain equivalent results, seal slides or keep humidified to minimize dehydration of mounting medium, store in dark at 2-8°C, and read within three days. Positive reactivity may range in fluorescence intensity from brilliant to weak. Grade the fluorescence reaction according to the following intensity scale: 4* (brilliant), 3* (bright), 2* (moderate) and 1* (weak).

Interpretation of Results

- Bright green fluorescent staining of the nuclei of the cells, in a clear pattern of fluorescence, denotes an ANA positive reaction. Samples that are positive at a 1:40 dilution should be titrated to the endpoint dilution. In a titration series, the highest

serum dilution demonstrating a 1⁺ reaction is interpreted as the endpoint. The observable patterns of nuclear fluorescence are:

Pattern	Characteristics	Associated Diseases
Homogeneous	Diffuse staining of the entire nucleus due to antibodies reactive with DNA-Nucleoprotein-Histone complexes (15).	SLE Rheumatoid Arthritis
Speckled	Specks of staining throughout the nucleus due to antibodies to Sm, RNP, SSA or SSB (16, 17).	Scleroderma Raynaud's Syndrome Sjögren's Syndrome Mixed Connective Tissue Disease
Nucleolar	Staining of the nucleolar structure due to antibodies against RNA-nucleoprotein complexes (10).	Scleroderma
Peripheral	Staining of the nuclear membranes due to antibodies against DNA (18, 19).	SLE
Anti-Centromere	Staining of the centromere regions of chromosomes due to reaction with a protein tightly bound to centromeric DNA (9, 10).	CREST Syndrome Variant of Progressive Systemic Sclerosis

- Absence of specific fluorescent pattern associated with the cell nucleus denotes an ANA negative reaction.

Significance of Interpretation

No discernible fluorescence of the cellular nuclei found at the screening dilution.	Test sample is negative for IgG ANA.
Specific positive fluorescence of the cellular nuclei found at the screening dilution or at higher dilutions.	Test sample is positive for IgG ANA.

Quality Control

- To ensure the test is working properly, use the positive and negative controls at least once for each day's testing. If the expected result for the controls is not met, the results are invalid and the test should be repeated.
- The type and age of the fluorescence microscope and hours of UV bulb usage can affect fluorescence intensity and titration endpoints to some degree. The ANA positive control furnished with this kit demonstrates a 3⁺ to 4⁺ intensity reaction. The vial has a listed titer to use as an additional check for the test system (see 1⁺ Dilution Notice). Use this as the calibrator for a 3⁺ to 4⁺ intensity reaction on your microscope.
- Use the ANA negative control furnished with this kit as the calibrator for a negative reaction on your equipment.
- Each day's test should include one PBS well in place of a test sample. This is a conjugate control to ensure the conjugate is not reacting with the cell substrate.
- Additional controls may be tested according to guidelines or requirements of local, state or federal regulations or accrediting organizations. For guidance on appropriate quality control practices, please refer to NCCLS document C24-A, Internal Quality Control Testing: Principles and Definitions.

1⁺ Dilution Notice

The positive control in this kit is ready to use and provides a 3⁺ to 4⁺ fluorescent intensity when tested. To obtain a 1⁺ fluorescent intensity make two-fold dilutions to the titer indicated on the vial included in this kit. Titer the positive control with the initial use of the kit.

The titer you obtain in testing may differ from the listed endpoint titer due to a number of technical reasons. It is best to test the titer indicated on the vial, as well as the two-fold dilution immediately preceding and following the listed titer. It is normal for results obtained for an endpoint (1⁺) titer to differ between laboratories due to factors affecting the intensity of fluorescence. These factors include:

- the power rating of the UV light source in the microscope
- the kind of light source
- the age of the lamp
- the length of the optical path of the microscope and the types of optical filters used
- the accuracy of dilution techniques and the dilution equipment

Limitations of the Procedure

- No single serological determination should be used as a sole criterion for detection. All available clinical and laboratory data should be included for diagnostic purposes.
- Patients with SLE or other autoimmune phenomena may show wide variations in ANA titers depending upon the clinical stage of the disease. An ANA titer is, therefore, not usually helpful in distinguishing disease states. High titers of ANA (1:1280 or above), however, may provide evidence of MCTD, SLE or scleroderma, especially if patterns of nucleolar, speckled or peripheral staining are observed.
- A number of therapeutic drugs have been reported to induce antinuclear antibodies in patients (e.g. procainamide, hydralazine, oral contraceptives, diphenylhydantoin, mephenytoin and some antibiotics). It is important to consider this information in determining evidence of disease.
- Low titers of antinuclear antibody may be found in the sera of normal individuals. Apparently healthy individuals may contain ANA in their sera. This percentage increases with age, particularly past the 6th decade of life (5, 22).
- If a sample is positive for more than one pattern, one of the patterns may partially or completely obscure the other. In this circumstance titration of the serum is required.
- SLE patients being treated with steroids may produce negative test results.
- The assay should be performed on serum. The use of other matrices has not been established.

Expected Values

The prevalence of ANA IgG antibodies in the normal population is expected to be zero. However, apparently healthy, asymptomatic individuals may contain ANA in their sera. This percentage increases with age, particularly past the sixth decade of life.

To demonstrate the prevalence of ANA in a normal population, sera from 99 normal individuals of various ages and genders, from various geographic locations were tested on the SCIMEDX anti-ANA IFA. Eighty-eight samples were negative for IgG antibodies to ANA (88.9%). Eleven samples were positive for IgG antibodies to ANA (11.1%).

Anti-ANA Prevalence in a Normal Population

Age Group	Total Samples	# Positive	Prevalence
12 - 19	20	1	5%
20 - 29	10	0	0%
30 - 39	16	2	12.5%
40 - 49	13	1	7.7%
50 - 59	10	0	0%

60 - 69	10	1	10.0%
>70	20	6	30.0%
Total	99	11	11.1%

To demonstrate the prevalence of ANA in a population of SLE patients' sera from 100 SLE patients were tested on the SCIMEDX anti-ANA IFA. Ninety-two samples were positive for IgG antibodies to ANA (92%). Eight samples were negative for IgG antibodies to ANA (8%).

Performance Characteristics

Relative Sensitivity and Specificity: The SCIMEDX ANA IFA kit was evaluated relative to a commercially available ANA IFA. The samples were frozen retrospective sera. One hundred sera were from patients diagnosed with SLE. Ninety-nine sera were from normals of various ages, gender and geographical areas. The overall agreement was 188/197 or 95%. The table below summarizes the data.

SCIMEDX ANA IFA				
Alternate IFA	ANA Status	Positive	Negative	Total
	Positive	100*	7	107
	Negative	2	88	90
	Total	102	95	197

*98/100 positive sera demonstrated the same pattern on the two assays. 2/100 samples showed different patterns.

Please be advised that 'relative' refers to the comparison of this assay's results to that of a similar assay. There was not an attempt to correlate the assay's results with disease presence or absence. No judgment can be made on the comparison assay's accuracy to predict disease.

Titer Agreement: Twenty-five positive sera were serially two-fold diluted and the endpoint titer was determined on the SCIMEDX ANA IFA and a commercially available ANA IFA. The endpoint titer results are as follows:

Identical Titer	12/25
± one, two-fold dilution	11/25
± two, two-fold dilutions	2/25

Reproducibility: Three positive sera with various titers (1:80, 1:640 and 1:2560) and one negative serum were serially diluted and assayed three times each on three different assays and the endpoint determined. All endpoint titers were within the specifications of ± one two-fold dilution. Refer to the following table.

Identical Titer	35/36
± one, two-fold dilution	1/36

Real Time Stability: Real time stability testing of kit components was carried out at 6 month intervals for a minimum of 24 months. The endpoint titers of the positive and negative controls were compared to the endpoint titers at release. Acceptance criteria are endpoint titers within one two-fold dilution of each other. These results were within specifications.

Real Time Stability

Slide Lot	Control	Endpoint Titer at Release	Endpoint Titer at 24 months
#1	Positive	1:1280	1:1280
	Negative	-	-
#2	Positive	1:1280	1:640
	Negative	-	-
#3	Positive	1:1280	1:640
	Negative	-	-

Literature Cited

1. Sinha, A.A., T. Lopez, and H.O. McDevitt. 1990. Autoimmune diseases: the failure of self tolerance. Science. 248:1380-1388.
2. Kotb, M. 1995. Infection and immunity: a story of the host, the pathogen, and the copathogen. Clin. Immunol. Immunopathol. 74:10-22.
3. Nygård, N.R. and B.D. Schwartz. 1993. Infection and autoimmunity. Adv. Int. Med. 38:337-359.
4. Ehrenstein, M. and D. Isenberg. 1991. Autoimmunity associated with infection: leprosy, acute rheumatic fever and Lyme disease. Curr. Opin. Immunol. 3:930-935.
5. Condemi, J.J., E.V. Barnett, E.C. Atwater, R.F. Jacob, E.S. Mongan, and J.H. Vaughan. 1965. The significance of antinuclear factors in rheumatoid arthritis. Arth. and Rheum. 8:1080-1093.
6. Coons, A.H. 1958. Fluorescent antibody methods, p. 399-422. In J.F. Danielli (ed.), General cytochemical methods, vol. 1. Academic Press, New York.
7. Friou, G.J. 1958. Identification of nuclear component of interaction of lupus erythematosus globulin and nuclei. J. Immunol. 80:476-481.
8. Friou G.J., S.C. Finch, and K.D. Detre. 1958. Interaction of nuclei and globulin from lupus erythematosus serum demonstrated with fluorescent antibody. J. Immunol. 80:324-329.
9. Tan, E.M. 1982. Autoantibodies to nuclear antigens (ANA): their immunobiology and medicine. Adv. Immunol. 33:167-240.
10. Tan, E.M., G.P. Rodnian, G. Ignacio, Y. Moroi, M.J. Fritzler, and C. Peebles. 1980. Diversity of anti-nuclear antibodies in progressive systemic sclerosis. Anti-centromere antibody and its relationship to CREST syndrome. Arth. and Rheum. 23:617-625.
11. Burnham, T.K. and P.W. Bush. 1974. Antinuclear antibodies. I. Patterns of nuclear immunofluorescence. J. Invest. Dermatol. 62:526-534.
12. Hall, A.P., W.A. Bardwil, T.B. Bayles, A.D. Mednis, and N. Galins. 1960. The relations between the antinuclear, rheumatoid and L.E.-cell factors in the systemic rheumatic diseases. N. Engl. J. Med. 263:769-775.
13. Pollack, V.E. 1964. Antinuclear antibodies in families of patients with systemic lupus erythematosus. N. Engl. J. Med. 271:165-171.
14. Raskin, J. 1964. Fluorescent antibody studies of certain dermatoses. Arch. Derm. 89:569-578.
15. Luciano, A. and N.F. Rothfield. 1973. Patterns of nuclear fluorescence and DNA-binding activity. Ann. Rheum. Dis. 32:337-341.
16. Friou, G.J. 1964. Immunofluorescence and antinuclear antibodies. Arth. and Rheum. 7:161-166.
17. Beck, J.W. 1962. Partial identification of the "speckled" nuclear antigen. Lancet 2: 241-242.
18. Tan, E.M. and H.G. Kunkel. 1966. Characteristics of a soluble nuclear antigen precipitating with sera of patients with systemic lupus erythematosus. J. Immunol. 96:464-471.
19. Beck, J.S. 1963. Autoantibodies to cell nuclei. Scot. Med. J. 8:373-388.
20. Casals, S.P., G.J. Friou, and L.L. Meyers. 1964. Significance of antibody to DNA in systemic lupus erythematosus. Arth. and Rheum. 7:379-390.
21. Anderson, J.R., K.G. Gray, J.S. Beck, W.W. Buchanan and A.J. McElhinney. 1962. Precipitating autoantibodies in the connective tissue diseases. Ann. Rheum. Dis. 21:360-369.
22. Whittingham, S., J. Irwin, and I.R. McKay. 1969. Autoantibodies in healthy subjects. Aust. Ann. Med. 18:130-134.

SCIMEDX CORPORATION
53 Richbenton Rd
Dover, NJ 07801 USA
Tel: 800.221.5598 Tel: 973.625.8822
Fax: 973.625.8796
www.scimedx.com

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