INDIRECT FLUORESCENCE ASSAY FOR HUMAN HERPESVIRUS 6 (HHV6) IgG ANTIBODY

EXPORT ONLY

I-HV601G	120 Tests
I-HV602G	40 Tests
I-HV603G	40 Tests

Intended Use

The SCIMEDX Indirect Fluorescence Assay (IFA) for Human Herpesvirus 6 Antigen (HHV6) IgG Antibody is intended for the qualitative and semi-quantitative detection of IgG (Immunoglobulin G) antibody to HHV6 in human serum or plasma. Detection of HHV6 IgG antibody in humans can be used as an aid in the diagnosis of primary infection or reactivation/reinfection with this virus.

Introduction and Summary of Test Procedures

Human herpesvirus 6 (HHV6) is a human lymphotropic virus, first isolated from immunosuppressed individuals and patients with lymphoproliferative disorders. It is a DNA virus, originally named human B-lymphotropic virus (HBLV), now classified as a member of the family of human herpesviruses (1, 2).

HHV6 is genetically and serologically distinct from the other known human herpesviruses (2-4). Serological studies, using the indirect fluorescence assay (IFA), have shown that antibody to the virus is often found in the general population. Primary infection usually occurs before the age of two. A high percentage of humans above the age of one are seropositive for anti-HHV6 antibody (4-6). Primary infection in adults is rare because of the high rate of infection in childhood. While primary infection in healthy adults usually results in mild illness, more serious complications can arise in immunocompromised patients or patients undergoing organ or bone marrow transplantation (4).

HHV6 has been identified as the causative agent of roseola infantum, also known as exanthem subitum. The virus has been isolated from children with exanthem subitum and the association further supported by serologic studies (7-9).

Roseola infantum is considered a benign non-seasonal disease of childhood with clinical indications of rash and fever. Infants are protected by maternal antibodies so that infection is most frequent between 6 and 18 months of age. Serological information is valuable in diagnosis since detection of antibody by the IFA technique is both sensitive and specific, without cross-reaction to other human herpesviruses (3, 5, 10, 11-13).

Principle of the Test

SCIMEDX fluorescent antibody assays use the indirect method of antibody detection and titer determination. Patient serum or plasma samples are applied to cultured cells containing inactivated viral antigens provided on paint delineated wells on glass microscope slides. During a 30 minute incubation, antibody specific for HHV6 antigens forms an antigen/antibody complex with the HHV6 antigens in the infected cells. In a brief washing step, nonspecific antibody and other unreacted serum proteins are eliminated. Fluorescein-conjugated goat anti-human IgG is then applied to the wells of the glass slide. The anti-IgG conjugate combines with human IgG, if present, during a 30 minute incubation. After a brief wash to remove unreacted conjugate, the slides are viewed by fluorescence microscopy. A positive antibody reaction is denoted by bright green fluorescence at the antigen sites.

Materials Furnished and Storage Conditions

HHV6 Antigen Slides: Slides of human lymphocytes infected with HHV6 on each glass well. The slides are ready for use after removal from protective pouch. Store at 2-8°C. Slides are stable at this storage condition until the expiration date stated on the pouch label.

HHV6 IgG Positive Control: Each vial contains 0.5 ml HHV6 IgG antibody positive human control. This component is a ready for use liquid. Store at 2-8°C. Liquid positive control is stable at this storage condition until the expiration date stated on the vial label.

HHV6 IgG Negative Control: Each vial contains 0.5 ml HHV6 IgG antibody negative human control. This component is a ready for use liquid. Store at 2-8°C. Liquid negative control is stable at this storage condition until the expiration date stated on the vial label.

Fluorescein Conjugate: Each vial contains 1.5 ml fluorescein conjugated goat (inactivated) antihuman IgG (heavy and light chain) with Evans Blue and Rhodamine counterstains. The fluorescein conjugate is a conjugation of affinity chromatography purified antihuman IgG with fluorescein isothiocyanate (FITC). Adding Evans Blue and Rhodamine counterstains to the conjugate masks nonspecific fluorescence of the tissue culture cells. This component is a ready for use liquid. Store at 2-8°C. Liquid conjugate is stable at this storage condition until the expiration date stated on the vial label.

Coverslip Mounting Media: Each vial contains 2.0 ml phosphate buffered glycerol with fade retardant. This component is a ready for use liquid. Storage temperature may range from refrigerated to room temperatures (2-30°C). Mounting media is stable at either storage condition until the expiration date stated on the vial label.

Phosphate Buffered Saline (PBS): Each aluminized sealed packet of powdered buffer makes one liter of 1X PBS. Storage temperatures may range from refrigerated to room temperature (2-30°C). Add the entire contents of a PBS packet to one liter of freshly prepared distilled or deionized water. Note: Addition of the salts while rapidly stirring the water will facilitate solubilization. Store PBS as a solution at 2-8°C.

Special Blotters: Absorbent blotters have pre-cut holes for use in drying the slide mask. Storage temperatures may range from refrigerated to room temperature (2-30°C).

General Precautions

IFA Test Kit: No US Standard of Potency.

- Store all kit components at their recommended or suggested temperatures. **Do not freeze.**
 - Do not use the components beyond the stated expiration date of each component.
 - SCIMEDX optimizes all of the active components in each lot of its IFA kits as a unit. Do not mix components from different lots or from different sources.
 - The controls and conjugate contain 0.095% sodium azide that, if allowed to accumulate, can form explosive compounds in lead and/or copper plumbing. Flush drains thoroughly if used to dispose of these materials.
 - Use separate pipette tips for each sample and reagent to avoid cross contamination.
 - Reagents should be inspected for evidence of bacterial or fungal contamination.
- Antigen slides:** All IFA antigen slides have fixed cells that do not contain any viable infectious agents. However, good laboratory practices (GLP) require careful slide handling and disposal as with any other potentially biohazardous laboratory material.
- Do not remove the slides from their protective pouch until ready for use.

- Do not reuse substrate slide.

Human Controls: The human controls in these kits have all been tested for hepatitis B surface antigen (HBsAg) and human immunodeficiency virus (HIV) antibody by FDA licensed methods and found to be non-reactive. However, no test system can ensure the absence of these agents. Handle all human serum components, including those received in your laboratory for testing, as potentially biohazardous.

Xn-Harmful Substance

Control and Conjugate Safety Precaution:

Concentration of sodium azide in these components is classified as Harmful and subject to the following risk phrase: "Harmful if swallowed and contact with acids liberates very toxic gas."

US

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		<ul style="list-style-type: none"> • Wash thoroughly after handling. • Wear eye/face protection
Signal Word	WARNING	<ul style="list-style-type: none"> • IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. • If eye irritation persists: Get medical advice/attention.
Hazard Statement	Causes serious eye irritation.	<ul style="list-style-type: none"> • Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

EU

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		<p>P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.</p>
Signal Word	WARNING	<p>P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.</p>
Hazard Statement	H319 Causes serious eye irritation.	<p>P337+P313 If irritation persists, get medical advice/attention.</p>

Specimen Collection, Storage and Limitations of Test Samples

- Separate aseptically collected serum or plasma from the red blood cells and store frozen (-10°C or lower) until ready for testing. Avoid repeated freezing and thawing.
- If desired, store fresh liquid serum or plasma samples at 2° to 8°C for up to one week without loss of antibody activity.
- Do not use excessively lipemic samples without delipidization.
- Do not use contaminated samples.
- Paired serum samples to demonstrate seroconversion or significant titer increase should be collected 7-14 days apart, stored, then tested simultaneously.

Additional Materials Required

- Test tubes, racks, pipettes, microtiter plates and safety pipetting devices for making sample dilutions
- 37°C incubator
- Moist chamber for incubating slides
- Slide holder rack and staining dish for washing slides
- Coverslips: 22 X 50 mm No. 1 thickness glass
- Fluorescence microscope: A fluorescence microscope equipped with the following was used to calibrate the controls and conjugate:
 - 10X eyepiece
 - 16X or 40X objectives
 - Epi-illuminator with 50W halogen lamp
 - FITC-excitation filter KP490
 - Yellow absorbing filter K530
 - Red suppression filter BG38

The fluorescein label has an excitation peak of 490 nm and an emission peak of 520 nm. Differences in endpoint reactivities and fluorescence intensities may be due to the type and condition of the fluorescence equipment used in your laboratory.

IFA Procedure

1. For qualitative IgG antibody determination, prepare a 1:20 screening dilution of each test sample in PBS. Prepare all dilutions in a minimum volume of 0.10 ml with PBS as the diluent.
2. For quantitative titration of sera, prepare two-fold serial dilutions of the serum sample in PBS, starting with a 1:20 dilution, and adding equal volumes of diluted serum or plasma and PBS for each consecutive dilution.
3. Remove slides from protective pouch and apply 1 drop (approximately 20 µl) of the diluted test sample(s) to each well. Add sufficient volume to completely cover each well, but cross-mixing of contents between wells should not occur. Note: Each day's test run requires one well each for positive control, negative control, and PBS (conjugate control).
4. Incubate the slides in a moist chamber for 30 minutes at 37°C.
5. Rinse the slides in a light stream of buffer. Avoid directing the stream at the wells.
6. Wash the slides for 10 minutes with a change of PBS solution after 5 minutes. Agitate the slides by moving the rack up and down in the buffer.
7. Blot the paint mask surrounding the test wells with the special blotters.
8. Apply one drop of the ready to use conjugate to each test well.
9. Repeat steps 4 (incubation), 5 (PBS rinse), 6 (10 minute PBS wash), and 7 (blot).
10. Apply the glycerol mounting media and 22 X 50 mm glass coverslip.
11. Observe the reactivity under fluorescence microscopy using 20-40X magnification. For best results, examine slides immediately after completion of the test. To obtain equivalent results, seal slides or keep humidified to minimize dehydration of mounting medium, store in dark at 2-8°C, and read within three days. Positive reactivity may range in fluorescence intensity from brilliant to weak. Grade the fluorescence reaction according to the following intensity scale: 4⁺ (brilliant), 3⁺ (bright), 2⁺ (moderate) and 1⁺ (weak).

Interpretation of Results

- Bright green fluorescent staining of the infected cells denotes a HHV6 IgG antibody positive reaction. The fluorescent staining pattern of HHV6 infected cells is variable. Depending on the cell's stage of infection, the fluorescent pattern can vary from a small portion of the infected cell fluorescing to the whole cell fluorescing. Fluorescence can

also range from granular to homogeneous. To provide an internal control, each well on the microscope slide contains both HHV6 infected and uninfected cells. Preparation of the slide in this manner is intentional. Uninfected cells, stained red by the counterstain, provide a contrasting background. Infectivity of the cells ranges from 20% - 60%. Titration of positive HHV6 IgG sera will provide quantitative information. In a titration series, the highest serum dilution demonstrating a 1⁺ reaction is interpreted as the endpoint.

- Absence of specific fluorescent staining of the infected cells denotes a HHV6 IgG antibody negative reaction.

Significance of Interpretation

1. No discernible fluorescence of the infected cells found at the screening dilution.	1. Test sample is HHV6 IgG antibody negative.
2. Specific positive fluorescence of the infected cells found at the screening dilution or at higher dilutions.	2. Test sample is HHV6 IgG antibody positive, indicating previous HHV6 infection. Seroconversion or a four-fold or greater rise in IgG antibody titer in paired serum samples indicates recent infection with HHV6.
3. Fluorescence found in both infected and uninfected cells.	3. Test sample is exhibiting a nonspecific reaction.

Note: Performance of a HHV6 IgM antibody specific test aids in the diagnosis of a current HHV6 infection.

Quality Control

- To ensure the test is working properly, use the positive and negative controls at least once for each day's testing.
- The type and age of the fluorescence microscope and hours of UV bulb usage can affect fluorescence intensity and titration endpoints to some degree. The HHV6 antibody positive control furnished with this kit demonstrates a 3⁺ to 4⁺ intensity reaction. The vial has a listed titer to use as an additional check for the test system (see 1⁺ Dilution Notice). Use this as the calibrator for a 3⁺ to 4⁺ intensity reaction on your microscope.
- Use the HHV6 antibody negative control furnished with this kit as the calibrator for a negative reaction on your equipment.
- Each day's test should include one PBS well in place of a test sample. This is a conjugate control to ensure the conjugate is not reacting with the cell substrate.

1⁺ Dilution Notice

The positive control in this kit is ready for use and provides a 3⁺ to 4⁺ fluorescent intensity when tested. To obtain a 1⁺ fluorescent intensity, make two-fold dilutions to the titer indicated on the vial included in this kit. Titer the positive control with the initial use of the kit. The titer you obtain in testing may differ from the listed dilution due to a number of technical reasons. It is best to test the titer indicated on the vial, as well as the two-fold dilution immediately preceding and following the listed titer. It is normal for results obtained for an endpoint (1⁺) titer to differ between laboratories due to factors affecting the intensity of fluorescence. These factors include:

- the power rating of the UV light source in the microscope
- the kind of light source
- the age of the lamp
- the length of the optical path of the microscope and the types of optical filters used
- the accuracy of dilution techniques and the dilution equipment

Limitations of the Procedure

- A serological test such as the IFA serves as an aid to detect viral infection, but its use should not be the sole criteria. The test results should be used in conjunction with information available from the patient, clinical evaluation and other available diagnostic procedures.
- A single positive result for HHV6 IgG antibody is significant only in that it indicates previous contact or infection with the virus. For epidemiological purposes a single result is useful. It should not be used, however, as an indication of current or recent infection with the virus.
- To determine current or recent infection, simultaneous testing of paired specimens of plasma or serum taken 7-14 days apart should be done. A four-fold or greater rise in titer between the first and second sample is indicative of a current or recent infection.
- Nonspecific positive reactions such as antinuclear antibody and/or anticytoplasmic antibody reactions can occur in samples from patients with certain autoimmune diseases. Both infected and uninfected cells will fluoresce, and this may obscure a positive HHV6 reaction. Therefore, observation of an autoimmune reaction cannot eliminate the possibility of HHV6 infection.

Expected Values

A high percentage of humans above the age of one are seropositive for anti-HHV6 antibody (4-6). Studies using the IFA method indicate prevalence rates of 80-95% in older children and adults in the general population. No significant differences in prevalence rates due to geographical location have been reported.

Performance Characteristics

Relative Sensitivity and Specificity: The SCIMEDX HHV6 IFA kit was evaluated relative to a commercially available HHV6 IFA. The samples were frozen retrospective sera. Fifty-seven sera were from donors with known HHV6 infection. Twenty-three sera were from normals of various ages, both genders and geographical areas. The overall agreement was 76/80 or 95.0%. Refer to the following table.

SCIMEDX HHV6 IFA

	HHV6 Status	Positive	Negative	Total
Alternate IFA	Positive	57	1	58
	Negative	3	19	22
	Total	60	20	80

Please be advised that 'relative' refers to the comparison of this assay's results to that of a similar assay. There was not an attempt to correlate the assay's results with disease presence or absence. No judgment can be made on the comparison assay's accuracy to predict disease.

Reproducibility: Ten positive sera with various titers (1:20-1:640) and five negative serum were serially diluted and assayed three times each and the endpoint determined. Forty-four of the endpoint titers were within the specifications of \pm one two-fold dilution. Refer to the following table.

Identical Titer	39/45
\pm one, two-fold dilution	5/45
\pm two, two-fold dilutions	1/45

Specificity: Thirty-four sera positive for antibody activity to diseases that have the potential for cross-reactivity to HHV6 were tested on this IFA kit. Because of the prevalence of HHV6 in patients worldwide, cross-reactivity was observed. The table below summarizes the data.

Cross Reactivity Data-SCIMEDX HHV6 IgG

Disease Type	Total Samples	Positive Result
Cytomegalovirus	18	11/18
Epstein-Barr virus	32	24/32
Herpes Simplex virus 1	26	17/26
Herpes Simplex virus 7	22	18/22
Varicella Zoster virus	10	5/10
Total	34	24/34

Real Time Stability: Real time stability testing of kit components was carried out at 6 month intervals for a minimum of 24 months. The endpoint titers of the positive and negative controls were compared to the endpoint titers at release. Acceptance criteria are endpoint titers within one two-fold dilution of each other. These results were within specifications. Refer to the table below.

Real Time Stability

Slide Lot	Control	Endpoint Titer at Release	Endpoint Titer at 24 months
#1	Positive	1:1280	1:1280
	Negative	-	-
#2	Positive	1:640	1:320
	Negative	-	-
#3	Positive	1:640	1:640
	Negative	-	-

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