INDIRECT FLUORESCENCE ASSAY FOR HUMAN HERPESVIRUS 7 (HHV7) IgG ANTIBODY

EXPORT ONLY

I-HV701G 120 Tests
I-HV702G 40 Tests

Intended Use

The SCIMEDX Indirect Fluorescence Assay (IFA) for Human Herpesvirus 7 (HHV-7) IgG Antibody is intended for the qualitative and semi-quantitative detection of IgG (Immunoglobulin G) antibody to HHV-7 in human serum or plasma. Detection of HHV-7 IgG antibody in humans can be used as an aid in the diagnosis of current infection with the virus and can provide serological evidence of prior infection.

Introduction and Summary of Test Procedures

Human herpesvirus 7 (HHV-7) is a ubiquitous human lymphotropic virus that was identified in 1990 from the T cells of a healthy individual (2). Classified as a member of the family of human herpesviruses, HHV-7 is closely related to HHV-6 and HCMV (1, 7). Like HHV-6, HHV-7 exhibits a tropism for T-cells and has been successfully propagated in a continuous human T-cell line, SupT1 (1). HHV-7 has been frequently isolated from saliva which indicates that the salivary gland may be the site of HHV-7 replication (2, 5, 8). Serological studies by IFA have shown that antibody to HHV-7 is common in the general population (1, 4, 8). HHV-7 primary infection does appear to occur later than the typical primary HHV-6 infection before the age of 2. The prevalence of HHV-7 antibody has been shown to increase gradually with age (1, 8). HHV-7 has been associated with exanthem subitum, an illness known to be caused by HHV-6 (3). No other association with disease has as yet been identified. Further studies are needed to clarify the characteristics of this virus, such as its etiologic role in disease and epidemiology.

The IFA has been shown to be an effective test for detecting HHV-7 IgG antibody. Seroconversion or a 4-fold or greater increase in IgG antibody in paired serum samples is indicative of current infection.

Principle of the Test

SCIMEDX fluorescent antibody assays use the indirect method of antibody detection and titer determination. Patient serum or plasma samples are applied to cultured cells containing inactivated viral antigens provided on paint delineated wells on glass microscope slides. During a 30 minute incubation, antibody specific for HHV-7 antigens forms an antigen/antibody complex with the HHV-7 antigens in the infected cells. In a brief washing step, nonspecific antibody and other unreacted serum proteins are eliminated. Fluorescein-conjugated goat anti-human IgG is then applied to the wells of the glass slide. The anti-IgG conjugate combines with human IgG, if present, during a 30 minute incubation. After a brief wash to remove unreacted conjugate, the slides are viewed by fluorescence microscopy. A positive antibody reaction is denoted by bright green fluorescence at the antigen sites.

Materials Furnished and Storage Conditions

HHV-7 Antigen Slides: Slides of human lymphocytes infected with HHV-7 on each glass well. The slides are ready for use after removal from protective pouch. Store at 2-8°C. Slides are stable at this storage condition until the expiration date stated on the pouch label.

HHV-7 IgG Positive Control: Each vial contains 0.5 ml HHV-7 IgG antibody positive human control. This component is a ready for use liquid. Store at 2-8°C. Liquid positive control is stable at this storage condition until the expiration date stated on the vial label.

HHV-7 IgG Negative Control: Each vial contains 0.5 ml HHV-7 IgG antibody negative human control. This component is a ready for use liquid. Store at 2-8°C. Liquid negative control is stable at this storage condition until the expiration date stated on the vial label.

Fluorescein Conjugate: Each vial contains 1.5 ml fluorescein conjugated goat (inactivated) antihuman IgG (heavy and light chain) with Evans Blue and Rhodamine counterstains. The fluorescein conjugate is a conjugation of affinity chromatography purified antihuman IgG with fluorescein isothiocyanate (FITC). Adding Evans Blue and Rhodamine counterstains to the conjugate masks nonspecific fluorescence of the tissue culture cells. This component is a ready for use liquid. Store at 2-8°C. Liquid conjugate is stable at this storage condition until the expiration date stated on the vial label.

Coverslip Mounting Media: Each vial contains 2.0 ml phosphate buffered glycerol with fade retardant. This component is a ready for use liquid. Storage temperature may range from refrigerated to room temperatures (2-30°C). Mounting media is stable at either storage condition until the expiration date stated on the vial label.

Phosphate Buffered Saline (PBS): Each aluminized sealed packet of powdered buffer makes one liter of 1X PBS. Storage temperatures may range from refrigerated to room temperature (2-30°C). Add the entire contents of a PBS packet to one liter of freshly prepared distilled or deionized water. Note: Addition of the salts while rapidly stirring the water will facilitate solubilization. Store PBS as a solution at 2-8°C.

Special Blotters: Absorbent blotters have pre-cut holes for use in drying the slide mask. Storage temperatures may range from refrigerated to room temperature (2-30°C).

General Precautions

IFA Test Kit: No US Standard of Potency.

- Store all kit components at their recommended or suggested temperatures. **Do not freeze.**
- Do not use the components beyond the stated expiration date of each component.
- SCIMEDX optimizes all of the active components in each lot of its IFA kits as a unit. Do not mix components from different lots or from different sources.
- The controls and conjugate contain 0.095% sodium azide that, if allowed to accumulate, can form explosive compounds in lead and/or copper plumbing. Flush drains thoroughly if used to dispose of these materials.
- Use separate pipette tips for each sample and reagent to avoid cross contamination.
- Reagents should be inspected for evidence of bacterial or fungal contamination.

Antigen slides: All IFA antigen slides have fixed cells that do not contain any viable infectious agents. However, good laboratory practices (GLP) require careful slide handling and disposal as with any other potentially biohazardous laboratory material.

- Do not remove the slides from their protective pouch until ready for use.
- Do not reuse substrate slide.

Human Controls: The human controls in these kits have all been tested for hepatitis B surface antigen (HBsAg) and human immunodeficiency virus (HIV) antibody by FDA licensed methods and found to be non-reactive. However, no test system can ensure the absence of

these agents. Handle all human serum components, including those received in your laboratory for testing, as potentially biohazardous.

Xn-Harmful Substance

Control and Conjugate Safety Precaution:

Concentration of sodium azide in these components is classified as Harmful and subject to the following risk phrase: "Harmful if swallowed and contact with acids liberates very toxic gas." US

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement
Pictogram		Prevention: <ul style="list-style-type: none"> • Wash thoroughly after handling. • Wear eye/face protection Response: <ul style="list-style-type: none"> • IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. • If eye irritation persists: Get medical advice/attention. Storage/Disposal: <ul style="list-style-type: none"> • Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.
Signal Word	WARNING	
Hazard Statement	Causes serious eye irritation.	

EU

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement
Pictogram		Prevention: <p>P264 Wash thoroughly after handling.</p> <p>P280 Wear protective gloves and clothing.</p> Response: <p>P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing.</p> <p>P337+P313 If irritation persists, get medical advice/attention.</p>
Signal Word	WARNING	
Hazard Statement	H319 Causes serious eye irritation.	

Specimen Collection, Storage and Limitations of Test Samples

- Separate aseptically collected serum or plasma from the red blood cells and store frozen (-10°C or lower) until ready for testing. Avoid repeated freezing and thawing.
- If desired, store fresh liquid serum or plasma samples at 2° to 8°C for up to one week without loss of antibody activity.
- Do not use excessively lipemic samples without delipidization.
- Do not use contaminated samples.
- Paired serum samples to demonstrate seroconversion or significant titer increase should be collected 7-14 days apart, stored, then tested simultaneously.

Additional Materials Required

- Test tubes, racks, pipettes, microtiter plates and safety pipetting devices for making sample dilutions
- 37°C incubator
- Moist chamber for incubating slides
- Slide holder rack and staining dish for washing slides
- Coverslips: 22 X 50 mm No. 1 thickness glass

Fluorescence microscope: A fluorescence microscope equipped with the following was used to calibrate the controls and conjugate:

- 10X eyepiece
- 16X or 40X objectives
- Epi-illuminator with 50W halogen lamp
- FITC-excitation filter KP490
- Yellow absorbing filter K530
- Red suppression filter BG38

The fluorescein label has an excitation peak of 490 nm and an emission peak of 520 nm. Differences in endpoint reactivities and fluorescence intensities may be due to the type and condition of the fluorescence equipment used in your laboratory.

IFA Procedure

1. For qualitative IgG antibody determination, prepare a 1:20 screening dilution of each test sample in PBS. Prepare all dilutions in a minimum volume of 0.10 ml with PBS as the diluent.
2. For quantitative titration of sera, prepare two-fold serial dilutions of the serum sample in PBS, starting with a 1:20 dilution, and adding equal volumes of diluted serum or plasma and PBS for each consecutive dilution.
3. Remove slides from protective pouch and apply 1 drop (approximately 20 µl) of the diluted test sample(s) to each well. Add sufficient volume to completely cover each well, but cross-mixing of contents between wells should not occur. Note: Each day's run requires one well each for positive control, negative control, and PBS (conjugate control).
4. Incubate the slides in a moist chamber for 30 minutes at 37°C.
5. Rinse the slides in a light stream of buffer. Avoid directing the stream at the wells.
6. Wash the slides for 10 minutes with a change of PBS solution after 5 minutes. Agitate the slides by moving the rack up and down in the buffer.
7. Blot the paint mask surrounding the test wells with the special blotters.
8. Apply one drop of the ready to use conjugate to each test well.
9. Repeat steps 4 (incubation), 5 (PBS rinse), 6 (10 minute PBS wash), and 7 (blot).
10. Apply the glycerol mounting media and 22 X 50 mm glass coverslip.
11. Observe the reactivity under fluorescence microscopy using 20-40X magnification. For best results, examine slides immediately after completion of the test. To obtain equivalent results, seal slides or keep humidified to minimize dehydration of mounting medium, store in dark at 2-8°C, and read within three days. Positive reactivity may range in fluorescence intensity from brilliant to weak. Grade the fluorescence reaction according to the following intensity scale: 4⁺ (brilliant), 3⁺ (bright), 2⁺ (moderate) and 1⁺ (weak).

Interpretation of Results

- Bright green fluorescent staining of the infected cells denotes a HHV-7 IgG antibody positive reaction. The fluorescent staining pattern of HHV-7 infected cells is variable. Depending on the cell's stage of infection, the fluorescent pattern can vary from a small portion of the infected cell fluorescing to the whole cell fluorescing. Fluorescence can also range from granular to homogeneous. To provide an internal control, each well on the microscope slide contains both HHV-7 infected and uninfected cells. Preparation of the slide in this manner is intentional. Uninfected cells, stained red by the counterstain, provide a contrasting background. Infectivity of the cells ranges from 20% - 60%. Titration of positive HHV-7 IgG sera will provide quantitative information. In a titration

series, the highest serum dilution demonstrating a 1⁺ reaction is interpreted as the endpoint.

- Absence of specific fluorescent staining of the infected cells denotes a HHV-7 IgG antibody negative reaction.

Significance of Interpretation

1. No discernible fluorescence of the infected cells found at the screening dilution.	1. Test sample is HHV-7 IgG antibody negative.
2. Specific positive fluorescence of the infected cells found at the screening dilution or at higher dilutions.	2. Test sample is HHV-7 IgG antibody positive, indicating previous HHV-7 infection. Seroconversion or a four-fold or greater rise in IgG antibody titer in paired serum samples indicates recent infection with HHV-7.
3. Fluorescence found in both infected and uninfected cells.	3. Test sample is exhibiting a nonspecific reaction.

Note: Performance of a HHV-7 IgM antibody specific test aids in the diagnosis of a current HHV-7 infection.

Quality Control

- To ensure the test is working properly, use the positive and negative controls at least once for each day's testing.
- The type and age of the fluorescence microscope and hours of UV bulb usage can affect fluorescence intensity and titration endpoints to some degree. The HHV-7 antibody positive control furnished with this kit demonstrates a 3⁺ to 4⁺ intensity reaction. The vial has a listed titer to use as an additional check for the test system (see 1⁺ Dilution Notice). Use this as the calibrator for a 3⁺ to 4⁺ intensity reaction on your microscope.
- Use the HHV-7 antibody negative control furnished with this kit as the calibrator for a negative reaction on your equipment.
- Each day's test should include one PBS well in place of a test sample. This is a conjugate control to ensure the conjugate is not reacting with the cell substrate.

1⁺ Dilution Notice

The positive control in this kit is ready to use and provides a 3⁺ to 4⁺ fluorescent intensity when tested. To obtain a 1⁺ fluorescent intensity make two-fold dilutions to the titer indicated on the vial included in this kit. Titer the positive control with the initial use of the kit.

The titer you obtain in testing may differ from the listed endpoint titer due to a number of technical reasons. It is best to test the titer indicated on the vial, as well as the two-fold dilution immediately preceding and following the listed titer. It is normal for results obtained for an endpoint (1⁺) titer to differ between laboratories due to factors affecting the intensity of fluorescence. These factors include:

- the power rating of the UV light source in the microscope
- the kind of light source
- the age of the lamp
- the length of the optical path of the microscope and the types of optical filters used
- the accuracy of dilution techniques and the dilution equipment

Limitations of the Procedure

- A serological test such as the IFA serves as an aid to detect viral infection, but its use should not be the sole criteria. The test results should be used in conjunction with information available from the patient, clinical evaluation and other available diagnostic procedures.
- A single positive result for HHV-7 IgG antibody is significant only in that it indicates previous contact or infection with the virus. For epidemiological purposes a single result is useful. It should not be used, however, as an indication of current or recent infection with the virus.

To determine current or recent infection, simultaneous testing of paired specimens of plasma or serum taken 7-14 days apart should be done. A four-fold or greater rise in titer between the first and second sample is indicative of a current or recent infection.

- Nonspecific positive reactions such as antinuclear antibody and/or anticytoplasmic antibody reactions can occur in samples from patients with certain autoimmune diseases. Both infected and uninfected cells will fluoresce, and this may obscure a positive HHV-7 reaction. Therefore, observation of an autoimmune reaction cannot eliminate the possibility of HHV-7 infection.

Performance Characteristics

Relative Sensitivity and Specificity: The SCIMEDX HHV7 IFA kit was evaluated relative to a commercially available HHV7 IFA. The samples were frozen retrospective sera. Twenty-eight sera were from normals of various ages, both genders and geographical areas. The overall agreement was 27/28 or 96.4%. Refer to the following table.

SCIMEDX HHV7 IFA				
Alternate IFA	HHV7 Status	Positive	Negative	Total
Positive	21	0	21	
Negative	1	6	7	
Total	22	6	28	

Please be advised that 'relative' refers to the comparison of this assay's results to that of a similar assay. There was not an attempt to correlate the assay's results with disease presence or absence. No judgment can be made on the comparison assay's accuracy to predict disease.

Reproducibility: Eight positive sera with various titers (1:20-1:160) and three negative serum were serially diluted and assayed three times each and the endpoint determined. All endpoint titers were within the specifications of \pm one two-fold dilution. Refer to the following table.

Identical Titer \pm one, two-fold dilution	28/33
	5/33

Specificity: The IFA test for HHV7 is a specific test for detecting antibodies to HHV7. Ten sera positive for antibody activity to diseases that have the potential for cross-reactivity to HHV7 were tested on this IFA kit. No cross-reactivity was observed for any of the samples. The table below summarizes the data.

Cross Reactivity Data-SCIMEDX HHV7 IgG		
Disease Type	Total Samples	Positive Result
Cytomegalovirus	7	0/7
Epstein-Barr virus	10	0/10

Herpes Simplex virus 1&2	9	0/9
Herpes Simplex virus 6	5	0/5
Herpes Simplex virus 8	3	0/3
Varicella Zoster virus	10	0/10
Total	10	0/10

Real Time Stability: Real time stability testing of kit components was carried out at 6 month intervals for a minimum of 24 months. The endpoint titers of the positive and negative controls were compared to the endpoint tiers at release. Acceptance criteria are endpoint titers within one two-fold dilution of each other. These results were within specifications. Refer to the table below.

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