INDIRECT FLUORESCENCE ASSAY FOR HUMAN HERPESVIRUS 8 (HHV8 LATENT) IgG ANTIBODY

EXPORT ONLY

I-HV801G 120 Tests
I-HV803G 40 Tests

Intended Use

The SCIMEDX Indirect Fluorescence Assay (IFA) for Human Herpesvirus 8 Antigen (HHV8 Latent) IgG Antibody is intended for the qualitative and semi-quantitative detection of IgG (Immunoglobulin G) antibody to HHV8 latent antigens in human serum or plasma. Detection of HHV8 IgG antibody in humans can be used as an aid in the diagnosis of primary infection or reactivation/reinfection with this virus, or can provide evidence of previous HHV8 infection.

Introduction and Summary of Test Procedures

Human herpesvirus 8 (HHV8), also known as Kaposi's Sarcoma Herpesvirus (KSHV), was identified in 1994 by Chang from KS lesions. Unique DNA sequences that closely resembled regions of the Herpesvirus saimiri and Epstein-Barr virus were discovered that exist in all forms of Kaposi's sarcoma worldwide, whether the KS is associated with human immunodeficiency virus (HIV) infection or not. Classical KS is a rare malignancy occurring in elderly Mediterranean men. In patients with acquired immunodeficiency syndrome (AIDS), KS is the most frequent neoplasm. KS has also been observed in transplant recipients and in certain African populations. In addition, HHV8 has also been associated with body cavity lymphomas, also called primary effusion lymphomas, multi-centric Castleman's disease, non-Hodgkin's lymphoma and multiple myeloma (1-6).

HHV8 is now classified as a gamma herpesvirus, in the genus Rhadinovirus, and resembles EBV in its tropism for B cells and ability to exist in a latent state. Although cell cultures established from body cavity lymphomas are often co-infected with HHV8 and EBV, some of the cultures are EBV negative. If the cell cultures are latently infected with HHV8, they can be chemically induced to produce actively replicating virus. Therefore, these HHV8, non-EBV infected lines, have made serological studies possible to detect antibodies to different viral proteins, including latent antigens, especially latency-associated nuclear antigens (LANA), and lytic antigens expressed by replicating virus (7-8).

Serological studies have shown that infection with HHV8 precedes the tumor development of clinical KS. Studies suggest that HHV8 is sexually transmitted and also indicate that HHV8 viral infection is not as common in the general population as the other herpesviruses. Different serological tests have shown varying prevalence rates for HHV8 antibodies in the normal adult population, from as low as 4% to 25%. More accurate determinations of prevalence rates remain to be established, with the use of a combination of tests measuring antibodies to both lytic and latent antigens and standardization of test systems. When both lytic and latent antibody tests are utilized, nearly all of KS patients and ~90% of HIV infected homosexual men show antibodies to HHV8 (9-14).

Among serological assays, the IFA has been shown to be a practical and accurate serological method for detection and semi-quantitation of IgG antibodies to HHV8. Because of the necessity for detection of antibodies to both lytic and latent antigens, SCIMEDX has two indirect fluorescence assays for determination of HHV8 antibodies, with two kinds of slides for differential detection of antibodies to lytic and latent antigens.

Principle of the Test

SCIMEDX fluorescent antibody assays use the indirect method of antibody detection and titer determination. Patient serum or plasma samples are applied to cultured cells containing inactivated viral antigens provided on paint delineated wells on glass microscope slides. During a 30 minute incubation, antibody specific for HHV8 antigens forms an antigen/antibody complex with the HHV8 antigens in the infected cells. In a brief washing step, nonspecific antibody and other unreacted serum proteins are eliminated. Fluorescein-conjugated goat anti-human IgG is then applied to the wells of the glass slide. The anti-IgG conjugate combines with human IgG, if present, during a 30 minute incubation. After a brief wash to remove unreacted conjugate, the slides are viewed by fluorescence microscopy. A positive antibody reaction is denoted by bright green fluorescence at the antigen sites.

Materials Furnished and Storage Conditions

HHV8-Latent Antigen Slides: Slides of human lymphocytes expressing HHV8 latent antigens on each glass well. The slides are ready for use after removal from protective pouch. Store at 2-8°C. Slides are stable at this storage condition until the expiration date stated on the pouch label.

HHV8 IgG Positive Control: Each vial contains 0.5 ml HHV8 IgG antibody positive human control. This component is a ready for use liquid. Store at 2-8°C. Liquid positive control is stable at this storage condition until the expiration date stated on the vial label.

HHV8 IgG Negative Control: Each vial contains 0.5 ml HHV8 IgG antibody negative human control. This component is a ready for use liquid. Store at 2-8°C. Liquid negative control is stable at this storage condition until the expiration date stated on the vial label.

Fluorescein Conjugate: Each vial contains 1.5 ml fluorescein conjugated goat (inactivated) antihuman IgG (heavy and light chain) with Evans Blue and Rhodamine counterstains. The fluorescein conjugate is a conjugation of affinity chromatography purified anti human IgG with fluorescein isothiocyanate (FITC). Adding Evans Blue and Rhodamine counterstains to the conjugate masks nonspecific fluorescence of the tissue culture cells. This component is a ready for use liquid. Store at 2-8°C. Liquid conjugate is stable at this storage condition until the expiration date stated on the vial label.

Coverslip Mounting Media: Each vial contains 2.0 ml phosphate buffered glycerol with fade retardant. This component is a ready for use liquid. Storage temperature may range from refrigerated to room temperatures (2-30°C). Mounting media is stable at either storage condition until the expiration date stated on the vial label.

Phosphate Buffered Saline (PBS): Each aluminized sealed packet of powdered buffer makes one liter of 1X PBS. Storage temperatures may range from refrigerated to room temperature (2-30°C). Add the entire contents of a PBS packet to one liter of freshly prepared distilled or deionized water. Note: Addition of the salts while rapidly stirring the water will facilitate solubilization. Store PBS as a solution at 2-8°C.

Special Blotters: Absorbent blotters have pre-cut holes for use in drying the slide mask. Storage temperatures may range from refrigerated to room temperature (2-30°C).

General Precautions

IFA Test Kit: No US Standard of Potency.

- Store all kit components at their recommended or suggested temperatures. **Do not freeze.**
- Do not use the components beyond the stated expiration date of each component.
- SCIMEDX optimizes all of the active components in each lot of its IFA kits as a unit. Do not mix components from different lots or from different sources.

- The controls and conjugate contain 0.095% sodium azide that, if allowed to accumulate, can form explosive compounds in lead and/or copper plumbing. Flush drains thoroughly if used to dispose of these materials.
- Use separate pipette tips for each sample and reagent to avoid cross contamination.
- Reagents should be inspected for evidence of bacterial or fungal contamination.

Antigen slides: All IFA antigen slides have fixed cells that do not contain any viable infectious agents. However, good laboratory practices (GLP) require careful slide handling and disposal as with any other potentially biohazardous laboratory material.

- Do not remove the slides from their protective pouch until ready for use.
- Do not reuse substrate slide.

Human Controls: The human controls in these kits have all been tested for hepatitis B surface antigen (HBsAg) and human immunodeficiency virus (HIV) antibody by FDA licensed methods and found to be non-reactive. However, no test system can ensure the absence of these agents. Handle all human serum components, including those received in your laboratory for testing, as potentially biohazardous.

Xn-Harmful Substance

Control and Conjugate Safety Precaution:

Concentration of sodium azide in these components is classified as Harmful and subject to the following risk phrase: "Harmful if swallowed and contact with acids liberates very toxic gas."

US

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		Response: <ul style="list-style-type: none"> • Wash thoroughly after handling. • Wear eye/face protection
Signal Word	WARNING	Storage/Disposal: <ul style="list-style-type: none"> • IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. • If eye irritation persists: Get medical advice/attention.
Hazard Statement	Causes serious eye irritation.	<ul style="list-style-type: none"> • Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

EU

Component	CS001 PBS Powder Packets O-GLB01Y Mounting Medium	Precautionary Statement Prevention:
Pictogram		Response: <ul style="list-style-type: none"> P264 Wash thoroughly after handling. P280 Wear protective gloves and clothing.
Signal Word	WARNING	Storage/Disposal: <ul style="list-style-type: none"> P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses if easy to do. Continue rinsing. P337+P313 If irritation persists, get medical advice/attention.
Hazard Statement	H319 Causes serious eye irritation.	

Specimen Collection, Storage and Limitations of Test Samples

- Separate aseptically collected serum or plasma from the red blood cells and store frozen (-10°C or lower) until ready for testing. Avoid repeated freezing and thawing.
- If desired, store fresh liquid serum or plasma samples at 2° to 8°C for up to one week without loss of antibody activity.
- Do not use excessively lipemic samples without delipidization.
- Do not use contaminated samples.
- Paired serum samples to demonstrate seroconversion or significant titer increase should be collected 7-14 days apart, stored, then tested simultaneously.

Additional Materials Required

- Test tubes, racks, pipettes, microtiter plates and safety pipetting devices for making sample dilutions
- 37°C incubator
- Moist chamber for incubating slides
- Slide holder rack and staining dish for washing slides
- Coverslips: 22 X 50 mm No. 1 thickness glass

Fluorescence microscope: A fluorescence microscope equipped with the following was used to calibrate the controls and conjugate:

- 10X eyepiece
- 16X or 40X objectives
- Epi-illuminator with 50W halogen lamp
- FITC-excitation filter KP490
- Yellow absorbing filter K530
- Red suppression filter BG38

The fluorescein label has an excitation peak of 490 nm and an emission peak of 520 nm. Differences in endpoint reactivities and fluorescence intensities may be due to the type and condition of the fluorescence equipment used in your laboratory.

IFA Procedure

1. For qualitative IgG antibody determination, prepare a 1:64 screening dilution of each test sample in PBS. Prepare all dilutions in a minimum volume of 0.10 ml with PBS as the diluent.
2. For quantitative titration of sera, prepare two-fold serial dilutions of the serum sample in PBS, starting with a 1:64 dilution, and adding equal volumes of diluted serum or plasma and PBS for each consecutive dilution.
3. Remove slides from protective pouch and apply 1 drop (approximately 20 µl) of the diluted test sample(s) to each well. Add sufficient volume to completely cover each well, but cross-mixing of contents between wells should not occur. Note: Each day's test run requires one well each for positive control, negative control, and PBS (conjugate control).
4. Incubate the slides in a moist chamber for 30 minutes at 37°C.
5. Rinse the slides in a light stream of buffer. Avoid directing the stream at the wells.
6. Wash the slides for 10 minutes with a change of PBS solution after 5 minutes. Agitate the slides by moving the rack up and down in the buffer.
7. Blot the paint mask surrounding the test wells with the special blotters.
8. Apply one drop of the ready to use conjugate to each test well.
9. Repeat steps 4 (incubation), 5 (PBS rinse), 6 (10 minute PBS wash), and 7 (blot).
10. Apply the glycerol mounting media and 22 X 50 mm glass coverslip.
11. Observe the reactivity under fluorescence microscopy using 20-40X magnification. For best results, examine slides immediately after completion of the test. To obtain equivalent results, seal slides or keep humidified to minimize dehydration of mounting medium, store in dark at 2-8°C, and read within three days. Positive reactivity may range in fluorescence intensity from brilliant to weak. Grade the

fluorescence reaction according to the following intensity scale: 4⁺ (brilliant), 3⁺ (bright), 2⁺ (moderate) and 1⁺ (weak).

Interpretation of Results

- Bright green fluorescent staining of the infected cells denotes a HHV8 IgG antibody positive reaction. The fluorescent staining pattern for antibody to latent HHV8 antigens is particulate fluorescence in the nucleus of the cells. To provide an internal control, each well on the microscope slide contains both HHV8 infected and uninfected cells. Preparation of the slide in this manner is intentional. Uninfected cells, stained red by the counterstain, provide a contrasting background. Infectivity of the cells ranges from 50-80%. Titration of positive HHV8 IgG sera will provide quantitative information. In a titration series, the highest serum dilution demonstrating a 1⁺ reaction is interpreted as the endpoint.
- Absence of specific fluorescent staining of the infected cells denotes a negative reaction for HHV8-Latent IgG antibody.

Significance of Interpretation

1. No discernible fluorescence of the infected cells found at the screening dilution.	1. Test sample is HHV8 IgG antibody negative.
2. Specific positive fluorescence of the infected cells found at the screening dilution or at higher dilutions.	2. Test sample is HHV8 IgG antibody positive, indicating previous HHV8 infection. Seroconversion or a four-fold or greater rise in IgG antibody titer in paired serum samples indicates recent infection with HHV8.
3. Fluorescence found in both infected and uninfected cells.	3. Test sample is exhibiting a nonspecific reaction.

Note: Performance of an HHV8 IgM antibody specific test aids in the diagnosis of a current HHV8 infection.

Quality Control

- To ensure the test is working properly, use the positive and negative controls at least once for each day's testing.
- The type and age of the fluorescence microscope and hours of UV bulb usage can affect fluorescence intensity and titration endpoints to some degree. The HHV8 antibody positive control furnished with this kit demonstrates a 2⁺ to 4⁺ intensity reaction. The vial has a listed titer to use as an additional check for the test system (see 1⁺ Dilution Notice). Use this as the calibrator for a 2⁺ to 4⁺ intensity reaction on your microscope.
- Use the HHV8 antibody negative control furnished with this kit as the calibrator for a negative reaction on your equipment.
- Each day's test should include one PBS well in place of a test sample. This is a conjugate control to ensure the conjugate is not reacting with the cell substrate.

1⁺ Dilution Notice

The positive control in this kit is ready to use and provides a 2⁺ to 4⁺ fluorescent intensity when tested. To obtain a 1⁺ fluorescent intensity, make two-fold dilutions to the titer indicated on the vial included in this kit. Titer the positive control with the initial use of the kit.

The titer you obtain in testing may differ from the listed dilution due to a number of technical reasons. It is best to test the titer indicated on the vial, as well as the two-fold dilution immediately preceding and following the listed titer. It is normal for results obtained for an endpoint (1⁺) titer to differ between laboratories due to factors affecting the intensity of fluorescence. These factors include:

- the power rating of the UV light source in the microscope
- the kind of light source
- the age of the lamp
- the length of the optical path of the microscope and the types of optical filters used
- the accuracy of dilution techniques and the dilution equipment

Limitations of the Procedure

- A serological test such as the IFA serves as an aid to detect viral infection, but its use should not be the sole criteria. The test results should be used in conjunction with information available from the patient, clinical evaluation and other available diagnostic procedures.
- A single positive result for HHV8 IgG antibody is significant only in that it indicates previous contact or infection with the virus. For epidemiological purposes a single result is useful. It should not be used, however, as an indication of current or recent infection with the virus.
- To determine current or recent infection, simultaneous testing of paired specimens of plasma or serum taken 7-14 days apart should be done. A four-fold or greater rise in titer between the first and second sample is indicative of a current or recent infection.
- Nonspecific positive reactions such as antinuclear antibody and/or anticytoplasmic antibody reactions can occur in samples from patients with certain autoimmune diseases. Both infected and uninfected cells will fluoresce, and this may obscure a positive HHV8 reaction. Therefore, observation of an autoimmune reaction cannot eliminate the possibility of HHV8 infection. In the test for latent antibody, the possibility of a positive ANA reaction being read as positive for latent antibody does exist. Comparison with the positive control, with its particulate reaction, or the testing of the sample with a test specific for ANA could be helpful to eliminate false positives.

Expected Values

When both lytic and latent antibody tests are utilized, nearly all of KS patients and ~90% of HIV infected homosexual men show antibodies to HHV8. Different serological tests, including the IFA, ELISA, and immunoblot assays, have shown varying prevalence rates for HHV8 antibodies in the normal adult population – from as low as 4% to 25%. Because of the different sensitivity and specificity of the various immunoassays, further testing with standardization of detection systems will be required to clarify prevalence rates. Seropositivity is not common before puberty.

Performance Characteristics

Relative Sensitivity and Specificity

This kit was compared with an ELISA kit commercially available for *in vitro* diagnostic use to determine the sensitivity and specificity between the two tests. The study tested 311 human serum samples that consisted of samples from patients with Kaposi's sarcoma, Human Herpes Virus 6 and multiple myelomas. Sera from normal blood donors were also included. The data is summarized in the following table.

SCIMEDX IFA

Samples	HHV8-Latent Status	Reactive	Non-reactive	Total
Kaposi's Sarcoma	Reactive	86	16	102
	Non-reactive	0	0	
Blood	Reactive	0	1	

Alternate ELISA	Donors	Non-reactive	2	191	194
Multiple Myelomas	Reactive	0	0	10	
	Non-reactive	0	10		
HHV6 Positive	Reactive	0	0	5	
	Non-reactive	0	5		
Total		88	223	311	

Serological Agreement: =292/311 = 93.9%

Reproducibility

Four seropositive samples with various titers (1:128 to 1:4096) and three seronegative samples were serially two-fold diluted and assayed four times each on two different days and the endpoint determined. Twenty-seven of the twenty-eight endpoint titers were within specification of plus or minus one two-fold dilution.

Identical titer	23/28
± one, two-fold dilution	4/28
± two, two-fold dilutions	1/28

HHV8 Specificity

The IFA test for HHV8-Latent is a specific test for detecting antibodies to HHV8. Thirty-four of thirty-five sera positive for IgG antibody to either cytomegalovirus (CMV), Epstein-Barr virus (EBV), human herpes virus 1 (HSV1), human herpes virus 6 (HHV6), or varicella zoster virus (VZV), or any combination of these viruses, tested negative for HHV8 IgG antibody with this IFA test. The table below summarizes the data.

Cross Reactivity Data-SCIMEDX HHV8 IgG

Disease Type	Total Samples	Positive Result
CMV	6	1/6
EBV	21	0/21
HSV1	1	0/1
HHV6	6	0/6
VZV	1	0/1
Total	35	1/35

Real Time Stability

Real time stability testing of kit components was carried out at 6 month intervals for a minimum of 24 months. The endpoint titers of the positive and negative controls were compared to the endpoint titers at release. Acceptance criteria are endpoint titers within one two-fold dilution of each other. These results were within specifications. Refer to the following table.

Real Time Stability

Slide Lot	Control	Endpoint Titer at Release	Endpoint Titer at 24 months
#1	Positive	1:512	1:256
	Negative	–	–
#2	Positive	1:512	1:512
	Negative	–	–
#3	Positive	1:512	1:256
	Negative	–	–

Literature Cited

- Chang, Y., E. Ceserman, M.S. Pessin, F. Lee, J. Culpepper, D.M. Knowles, and P.S. Moore. 1994. Identification of herpesvirus-like DNA sequences in AIDS-associated Kaposi's sarcoma. *Science* **265**: 1865-1869.
- Levy, J.A. 1997. Three new human herpesviruses (HHV6, 7 and 8). *Lancet* **349**: 558-562.
- Beral, V. 1996. Epidemiology of Kaposi's sarcoma, in *Cancer Surveys*, vol. 10: Cancer, HIV, and AIDS (Beral, V., Jaffe, H.W., Weiss, R., eds) pp5-22. Cold Spring Harbor Press, Cold Spring Harbor, New York.
- Ganem, D. 1996. Human herpesvirus 8 and the biology of Kaposi's sarcoma. *Virology* **7**: 325-332.
- Luppi, M., and G. Torrelli. 1996. The new lymphotropic herpesviruses (HHV-6, HHV-7, HHV8) and Hepatitis C Virus (HCV) in human lymphoproliferative diseases: an overview. *Haematologica* **81**: 265-281.
- Kim, S., S.H. Kang, J. Moon, M.S. Kim, Y.S. Kim, K. Park. 1998. Kaposi's sarcoma in renal transplant recipients. *Transplantation proceedings* **30**: 3165.
- Miller, G., M. Rigby, L. Heston, E. Grogan, R. Sun, C. Metroka, J. Levy, S-J. Gao, Y. Chang, P. Moore. 1995. Antibodies to butyrate-inducible antigens of Kaposi's sarcoma-associated herpesvirus in patients with HIV-1 infection. *New Engl J Med* **334**:1291-1297.
- Lennette, E.T., D.J. Blackburn, and J.A. Levy. 1996. Antibodies to human herpesvirus type 8 in the general population and in Kaposi's sarcoma patients. *Lancet* **348**: 858866.
- Blackburn, D.J., J. Ambroziak, E. Lennette, M. Adams, B. Ramachandran, and J.A. Levy. 1997. Infectious human herpesvirus 8 in a healthy North American blood donor. *Lancet* **349**: 609-611.
- Simpson, G.R., T.F. Schultz, D. Whitley, P.M. Cook, C. Boshoff, L. Rainbow, M.R. Howard, S.Gao, R.A. Bohenzky, P. Simmonds, C. Lee, A. de Ruiter, A. Hatzakis, R.S. Tedder, I.V.D. Weller, R.A. Weiss, and P.S. Moore. 1996. Prevalence of Kaposi's sarcoma associated herpesvirus infection measured by antibodies to recombinant capsid protein and latent immunofluorescence antigen. *Lancet* **349**: 1133-1138.
- Kedes, D.H., D. Ganem, N. Ameli, P. Bacchetti, and R. Greenblatt. 1997. The prevalence of serum antibody to human herpesvirus 8 (Kaposi's Sarcoma-Associated Herpesvirus) among HIV-seropositive and high-risk HIV-seronegative women. *JAMA* **277**: 478-481.
- Kedes, D.H., E. Operksalski, M. Busch, K. Kohn, J. Flood, D. Ganem. 1996. The seroepidemiology of human herpesvirus 8 (Kaposi's sarcoma-associated herpesvirus): distribution of infection in KS risk groups and evidence for sexual transmission. *Nature Med.* **2**: 918-921.
- Martin, J.N., D.Ee. Ganem, D.H. Osmond, K.A. Page-Shafer, D. Macrae, and D.H. Kedes. 1998. Sexual transmission and the natural history of human herpesvirus 8 infection. *New Engl. J Med.* **338**: 948-954.
- Chatlynne, L.G., W. Lapps, M. Handy, Y.Q. Huang, R. Masood, A.S. Hamilton, J.W. Said, H.P. Koeflffer, M.H. Kaplan, A. Friedman-Kien, P.S. Gill, J.E. Whitman, and D.V. Abashli. 1998. Detection and titration of human herpesvirus-8 specific antibodies in sera from blood donors, acquired immunodeficiency syndrome patients, and Kaposi's sarcoma patients using a whole virus enzyme-linked immunosorbent assay. *Blood* **92**: 53-58.

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
Literaturnachweis

1. Chang, Y., E. Cesarman, M.S. Pessin, F. Lee, J. Culpepper, D.M. Knowles, and P.S. Moore. 1994. Identification of herpesvirus-like DNA sequences in AIDS-associated Kaposi's sarcoma. *Science* **265**: 1865–1869.
2. Levy, J.A. 1997. Three new human herpesviruses (HHV6, 7 and 8). *Lancet* **349**: 558–562.
3. Beral, V. 1996. Epidemiology of Kaposi's sarcoma, in Cancer Surveys, vol. 10: Cancer, HIV, and AIDS (Beral,V., Jaffe, H.W., Weiss, R.,eds) pp5–22. Cold Spring Harbor Press, Cold Spring Harbor, New York.
4. Ganem, D. 1996. Human herpesvirus 8 and the biology of Kaposi's sarcoma. *Virology* **7**: 325–332.
5. Luppi, M., and G. Torelli. 1996. The new lymphotropic herpesviruses (HHV-6, HHV-7, HHV8) and Hepatitis C Virus (HCV) in human lymphoproliferative diseases: an overview. *Haematologica* **81**: 265–281.
6. Kim, S., S.H. Kang, J. Moon, M.S. Kim, Y.S. Kim, K. Park. 1998. Kaposi's sarcoma in renal transplant recipients. *Transplantation proceedings* **30**: 3165.
7. Miller, G. M. Rigsby, L. Heston, E. Grogan, R. Sun, C. Metroka, J. Levy, S.-J. Gao, Y. Chang, P. Moore. 1995. Antibodies to butyrate-inducible antigens of Kaposi's sarcoma-associated herpesvirus in patients with HIV-1 infection. *New Engl J Med* **334**:1291–1297.
8. Lennette, E.T., D.J. Blackburn, and J.A. Levy. 1996. Antibodies to human herpesvirus type 8 in the general population and in Kaposi's sarcoma patients. *Lancet* **348**: 858866.
9. Blackburn, D.J., J. Ambroziak, E. Lennette, M. Adams, B. Ramachandran, and J.A. Levy. 1997. Infectious human herpesvirus 8 in a healthy North American blood donor. *Lancet* **349**: 609–611.
10. Simpson, G.R., T.F. Schultz, D. Whity, P.M. Cook, C. Boshoff, L. Rainbow, M.R. Howard, S.Gao, R.A. Bohenzl, P. Simmonds, C. Lee, A. de Ruiter, A. Hatzakis, R.S. Tedder, I.V.D Weller, R.A Weiss, and P.S. Moore. 1996. Prevalence of Kaposi's sarcoma associated herpesvirus infection measured by antibodies to recombinant capsid protein and latent immunofluorescence antigen. *Lancet* **349**: 1133–38.
11. Kedes, D.H., D. Ganem, N. Armeli, P. Bacchetti, and R. Greenblatt. 1997. The prevalence of serum antibody to human herpesvirus 8 (Kaposi Sarcoma-Associated Herpesvirus) among HIV-seropositive and high-risk HIV-seronegative women. *JAMA* **277**: 478–481.
12. Kedes, D.H., E. Operksalski, M. Busch, K. Kohn, J. Flood, D. Ganem. 1996. The seroepidemiology of human herpesvirus 8 (Kaposi's sarcoma-associated herpesvirus): distribution of infection in KS risk groups and evidence for sexual transmission. *Nature Med.* **2**: 918–921.
13. Martin, J.N., D.Ee. Ganem, D.H. Osmond, K.A. Page-Shafer, D. Macrae, and D.H. Kedes. 1998. Sexual transmission and the natural history of human herpesvirus 8 infection. *New Engl. J Med.* **338**: 948–54.
14. Chatlynne, L.G., W. Lapps, M. Handy, Y.Q. Huang, R. Masood, A.S. Hamilton, J.W. Said, H.P. Koeffler, M.H. Kaplan, A. Friedman-Kien, P.S. Gill, J.E. Whitman, and D.V. Ablashi. 1998. Detection and titration of human herpesvirus-8 specific antibodies in sera from blood donors, acquired immunodeficiency syndrome patients, and Kaposi's sarcoma patients using a whole virus enzyme-linked immunosorbent assay. *Blood* **92**: 53–58.

 SCIMEDIX CORPORATION
53 Richboynton Rd
Dover, NJ 07801 USA
Tel: 800.221.5588 Tel: 973.625.8822
Fax: 973.625.8796
www.scimedix.com

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
7. Miller, G. M. Rigsby, L. Heston, E. Grogan, R. Sun, C. Metroka, J. Levy, S-J. Gao, Y. Chang, P. Moore. 1995. Antibodies to butyrate-inducible antigens of Kaposi's sarcoma-associated herpesvirus in patients with HIV-1 infection. *New Engl J Med* 334:1291-1297.
8. Lennette, E.T., D.J. Blackbourn, and J.A. Levy. 1996. Antibodies to human herpesvirus type 8 in the general population and in Kaposi's sarcoma patients. *Lancet* 348: 858866.
9. Blackbourn, D.J., J. Ambroziak, E. Lennette, M. Adams, B. Ramachandran, and J.A. Levy. 1997. Infectious human herpesvirus 8 in a healthy North American blood donor. *Lancet* 349: 609-611.
10. Simpson, G.R., T.F. Schultz, D. Whitby, P.M. Cook, C. Boshoff, L. Rainbow, M.R. Howard, S.Gao, R.A. Bohenzky, P. Simmonds, C. Lee, A. de Ruiter, A. Hatzakis, R.S. Tedder, I.V.D. Weller, R.A Weiss, and P.S. Moore. 1996. Prevalence of Kaposi's sarcoma associated herpesvirus infection measured by antibodies to recombinant capsid protein and latent immunofluorescence antigen. *Lancet* 349: 1133-38.
11. Kedes, D.H., D. Ganem, N. Ameli, P. Baccetti, and R. Greenblatt. 1997. The prevalence of serum antibody to human herpesvirus 8 (Kaposi Sarcoma-Associated Herpesvirus) among HIV-seropositive and high-risk HIV-seronegative women. *JAMA* 277: 478-481.
12. Kedes, D.H., E. Operksalski, M. Busch, K. Kohn, J. Flood, D. Ganem. 1996. The seroepidemiology of human herpesvirus 8 (Kaposi's sarcoma-associated herpesvirus): distribution of infection in KS risk groups and evidence for sexual transmission. *Nature Med.* 2: 918-921.
13. Martin, J.N., D.Ee. Ganem, D.H. Osmond, K.A. Page-Shafer, D. Macrae, and D.H. Kedes. 1998. Sexual transmission and the natural history of human herpesvirus 8 infection. *New Engl. J Med.* 338: 948-54.
14. Chatlynne, L.G., W. Lapps, M. Handy, Y.Q. Huang, R. Masood, A.S. Hamilton, J.W. Said, H.P. Koeffler, M.H. Kaplan, A. Friedman-Kien, P.S. Gill, J.E. Whitman, and D.V. Ablashi. 1998. Detection and titration of human herpesvirus-8 specific antibodies in sera from blood donors, acquired immunodeficiency syndrome patients, and Kaposi's sarcoma patients using a whole virus enzyme-linked immunosorbent assay. *Blood* 92: 53-58.

 SCIMEDIX CORPORATION
53 Richbenton Rd
Dover, NJ 07801 USA
Tel: 800.221.5598 Tel: 973.625.8822
Fax: 973.625.8796
www.scimedix.com

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- Weller, R.A Weiss, and P.S. Moore.** 1996. Prevalence of Kaposi's sarcoma associated herpesvirus infection measured by antibodies to recombinant capsid protein and latent immunofluorescence antigen. *Lancet* **349:** 1133-38.
11. **Kedes, D.H., D. Ganem, N. Ameli, P. Bacetti, and R. Greenblatt.** 1997. The prevalence of serum antibody to human herpesvirus 8 (Kaposi Sarcoma-Associated Herpesvirus) among HIV-seropositive and high-risk HIV-seronegative women. *JAMA* **277:** 478-481.
 12. **Kedes, D.H., E. Operalski, M. Busch, K. Kohn, J. Flood, D. Ganem.** 1996. The seroepidemiology of human herpesvirus 8 (Kaposi's sarcoma-associated herpesvirus); distribution of infection in KS risk groups and evidence for sexual transmission. *Nature Med.* **2:** 918-921.
 13. **Martin, J.N., D.E. Ganem, D.H. Osmond, K.A. Page-Shafer, D. Macrae, and D.H. Kedes.** 1998. Sexual transmission and the natural history of human herpesvirus 8 infection. *New Engl. J Med.* **338:** 948-54.
 14. **Chatlynne, L.G., W. Lapps, M. Handy, Y.Q. Huang, R. Masood, A.S. Hamilton, J.W. Said, H.P. Koeffler, M.H. Kaplan, A. Friedman-Kien, P.S. Gill, J.E. Whitman, and D.V. Ablashi.** 1998. Detection and titration of human herpesvirus-8 specific antibodies in sera from blood donors, acquired immunodeficiency syndrome patients, and Kaposi's sarcoma patients using a whole virus enzyme-linked immunosorbent assay. *Blood* **92:** 53-58.

 SCIMEDIX CORPORATION
53 Richboynton Rd
Dover, NJ 07801 USA
Tel: 800.221.5598 Tel: 973.625.8822
Fax: 973.625.8795
www.scimedix.com

EC	REP	Medimark Europe 11, rue Zola – BP 2332 F-38033 Grenoble Cedex 2 - France
----	-----	---

	Manufactured by Prodotto da Fabricado por Fabriqué par hergestellt von
REF	Catalog number Número de catalogo Número de Catálogo Numéro de catalogue Katalognummer
LOT	Lot Lotto Lote Lot Charge
EC REP	EC Authorized Representative Rappresentante Autorizzato CE Representante Autorizado CE CE Représentant autorisé EG autorisierter Bevollmächtigter
CE	EC Declaration of Conformity Dichiarazione di Conformità CE Declaración de Conformidad CE CE Déclaration de Conformité EG Konformitätserklärung
	Number of tests Número di test Número de determinaciones Nombre de tests Anzahl der Teste
	See instructions for use Vedere le istruzioni per l'uso Consultar la instrucciones de uso Voir instructions Gebrauchsweisung beachten
	Expiration date Data di scadenza Caducidad Date d'Expiration Haltbarkeitsdatum
	Store at 2-8 °C / 35-46 °F Conservare a 2-8°C/35-46 F Almacenar a 2-8°C / 35-46°F Conserver à 2-8°C Bei 2-8°C / 35-46°F lagern
	Caution Attenzione Precaución Précautions Achtung
	Xn-Harmful Substance. Refer to Material Safety Data Sheet Substance nocive (Xn). Voir la fiche sur la sécurité des équipements. Xn – gefährlicher Stoff Siehe Sicherheitsdatenblatt. Xn-Sostanza nociva. Fare riferimento alla Scheda di sicurezza. Xn - Sustancia dañina. Consultar la Hoja de Datos de Seguridad del Material.
	Potential biological risk Potenziale rischio biologico Riesgo potencial biológico Biohazard Potentielle biologische Gefährdung
RFU	Ready for use Pronto all'uso Listo para su uso Prêt à l'emploi Gebrauchsfertig

IVD	For in vitro diagnostic use Per uso diagnostico <i>in vitro</i> Para uso solo <i>in vitro</i> Usage <i>in vitro</i> Für in-vitro diagnostische Verwendung
RUO	For research use only Solo per ricerca Para uso solo en investigación Pour recherche Nur für Forschungszwecke
IUO	For investigational use only Solo per uso investigativo Para uso solo en investigación Pour investigation Nur für Forschungszwecke
PBS	Phosphate Buffered Saline Tampone salino fosfato Fosfato Salino Tamponado Tampon phosphate saline PBS
MS	Mounting Solution, Buffered Glycerol Solution de montage, glycérol tamponné Eindeckmittel, gepuffertes Glycerol Soluzione di montaggio, glicerolo tamponato Solución de montaje, glicerol tamponado
AG SLD	Antigen Slide Lames d'antigène Antigen-Objekträger: Vetrino dell'antigene Portaobjetos de antígeno
FITC	Anti-Human Conjugate with counterstain Conjugué anti-humain avec coloration de contraste Antihuman-Konjugat mit Gegenfärbung Conjugado anti-umano con colorante de contraste Conjugado anti-humano con contratinación
NEG	Negative Control Controllo Negativo Control Negativo Contrôle negative Negative Kontrolle
POS	Positive Control Controllo Positivo Control Positivo Contrôle positif Positive Kontrolle
CONJ CNS	Counterstain Colorante di contrasto Contraste Contrecolorant Gegenfärbung
CS	Coverslip Coprioggetto Cubre portaobjetos Lamelles couvre-objet Deckgläshen
BLT	Blotters Buvards Löschkopierpapier Tamponi di carta assorbente Absorbentes
DIL	Sample Diluent Diluente del campione Diluyente de muestra Tampon de dilution Probenverdünnungslösung